




Installing square optical elements, right-handed system
 Accepts a height of 25.4 mm (1 inch)
 Minimum thickness is 5 mm [0.2 inch]
 The maximum thickness is 6 mm [0.24 inch].




6.00
[0.24 in]
Maximum lens thickness




2-M4 cylindrical head screw
or M6 set screw
is used for fixing the spectacle frame.

Cross section A-A

3-M6X0.25 adjusting screw,
 Adjustable angle range ± 4 ,
 0.39 degrees per turn,
 With knob,
 2.0mm hexagon wrench hole,



DRAWING PROJECTION 			EMR25.4—3AR	Zolix				
	NAME	DATE		Rectangular adjusting frame	MATERIAL	WEIGHT	SCALE	REV
DRAWN	NIC	10th/June/2025			FOLLOW THE PARTS	0.1	1:1.2	V1
APPROVAL	DICK	10th/June/2025						
THE DRAWINGS ARE ONLY FOR TYPE SELECTION INFORMATION				shop.zolix.com.cn				