




Installing square optical elements, right-handed system
 Accepts a height of 50.8 mm (2 inches)
 Minimum thickness is 5 mm [0.2 inch]
 The maximum thickness is 7.5 mm [0.3 inch].




2-M4 cylindrical head screw
or M6 set screw
is used for fixing the spectacle frame.



Cross section A-A

2-M6X0.25 adjusting screw,
Adjustable angle range $\pm 4^\circ$,
0.25 degrees per lap,
With knob,
2.0mm hexagon wrench hole,



DRAWING PROJECTION 			EMR50.8-2AR		Zolix			
	NAME	DATE	Rectangular adjusting frame		MATERIAL	WEIGHT	SCALE	REV
DRAWN	NIC	10th/June/2025						
APPROVAL	DICK	10th/June/2025						
THE DRAWINGS ARE ONLY FOR TYPE SELECTION INFORMATION					FOLLOW THE PARTS	0.13	1:1.8	V1
					shop.zolix.com.cn			