




Installing square optical elements, right-handed system
 Accepts a height of 50.8 mm (2 inches)
 Minimum thickness is 5 mm [0.2 inch]
 The maximum thickness is 7.5 mm [0.3 inch].




2-M4 cylindrical head screw
or M6 set screw
is used for fixing the spectacle frame.



Cross section A-A

3-M6X0.25 adjusting screw,
Adjustable angle range $\pm 4^\circ$,
0.25 degrees per lap,
With knob,
2.0mm hexagon wrench hole,



DRAWING PROJECTION 			EMR50.8—3AR	Zolix				
	NAME	DATE		Rectangular adjusting frame	MATERIAL	WEIGHT	SCALE	REV
DRAWN	NIC	10th/June/2025			FOLLOW THE PARTS	0.14	1:1.6	V1
APPROVAL	DICK	10th/June/2025	shop.zolix.com.cn					
THE DRAWINGS ARE ONLY FOR TYPE SELECTION INFORMATION								